TO SHOW THAT a ∞ F

Apparatus

Air-track with one vehicle, pulley and blower, two photogates, two retort stands, dual timer, metre-stick, black card, set of slotted weights (1 N total).



Procedure

- 1. Set up the apparatus as in the diagram. Make sure the card cuts both light beams as it passes along the track.
- 2. Level the air track.
- 3. Set the weights F at 1 N. With the card at one end of the track start the blower and release the card from rest.
- 4. Note the times t_1 and t_2 .
- 5. Remove one 0.1 N disc from the slotted weight, store this on the vehicle, and repeat.
- 6. Continue for values of F from 1.0 N to 0.1 N.
- 7. Use a metre-stick to measure the length of the card *l* and the separation of the photogate beams *s*.
- 8. Record results as shown.
- 9. Draw a graph of $a/m \text{ s}^{-2}$ against F/N.

Results

l = m.

 $s = \dots m$.

Initial velocity $u = \frac{l}{t_1}$

Final velocity $v = \frac{l}{t_2}$

Acceleration $a = \frac{v^2 - u^2}{2s}$

F/N	t_1/s	t_2/s	$u/m s^{-1}$	$v/m s^{-1}$	$a/m s^{-2}$
1.0					
0.9					
0.8					
0.7					
0.6					
0.5					
0.4					
0.3					
0.2					
0.1					

Conclusion

A straight line through the origin shows that, for a constant mass, the acceleration is proportional to the applied force.

Notes

The total accelerating mass must be kept constant; hence the need to transfer the masses. Block the ten pairs of air holes nearest the buffer/pulley end of the track with cellotape. This part of the track will now act as a brake on the vehicle.

Occasionally check the air holes on the linear air-track with a pin, to clear any blockages due to grit or dust.

This experiment may be performed using a trolley on a friction-compensated ramp.