


# INVESTIGATION OF THE HETEROGENEOUS REACTION OF AMMONIA ON DELIQUESCENT SEASALT USING FTIR SPECTROSCOPY.

Arnaud Allanic, Colette Noonan, John Sodeau.

Centre for Research into Atmospheric Chemistry, University College Cork, and Environmental Research Institute, Cork, Ireland.  
Crac@ucc.ie



## Introduction



- Ammonia has mainly agricultural sources (live stock, fertiliser). It is expected to represent a significant atmospheric pollutant in Ireland (mixing ratio >50ppb)
- Adsorbed ammonia enhances the rate of oxidation of SO<sub>2</sub> to yield stable NH<sub>4</sub>HSO<sub>4</sub> or (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> aerosols.

- Hygroscopic behaviour of (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> displays a much stronger hysteresis than that of H<sub>2</sub>SO<sub>4</sub>.



- Uptake of ammonia into acidic aerosol can affect:
  - growth pattern and ability to act as cloud nuclei
  - alter optical properties (light absorption and scattering by particulate matter).

- Ammonia uptake rate on deliquescent aerosol is controlled by diffusion towards counter ions in the bulk as well as the kinetics of surface regeneration.

- Spectroscopic technique can provide a more complete picture of the mechanism, as it probes bulk phase processes.



## Experimental Setup



The aerosol suspension in nitrogen is generated by a TSI 3076 atomizer from a 1% wt seasalt solution. A diffusion dryer allows for complete re-crystallization of the salt dissolved in the droplets. The aerosol fraction is monitored downstream by a particle sizer (SMPS TSI 3081). The relative humidity (RH) in the reactor can be tuned between 1% and 90%.


### Laminar flow conditions

(Reynolds Number  $Re = 60$  at 5L/min.)

$$t = \frac{z}{f} \quad f: \text{flow velocity (cm s}^{-1}\text{)}$$

### First-order reaction rate k:


$$f \left( \frac{d[\text{NH}_3]}{dz} \right) = -k[\text{NH}_3]$$




The flow-tube is made of glass and has an inner diameter of 10cm and a maximal reactive length of 80cm. This reactive length - i.e. the fraction of the reactor located between the tip of the injector and the sampling inlet - can be tuned up stepwise to this maximal value using a movable injector (6mm outer diameter).

It operates at atmospheric pressure and its temperature can be regulated. It is fitted with spectroscopic windows (BaF<sub>2</sub>) to allow for FTIR monitoring of both gas- and condensed phase (Digilab FTS 3000 coupled with external MCT detector). The geometry of the reactor ensures that laminar flow conditions are always achieved, so that reaction time  $t$  is directly proportional to reaction length  $z$ .

## NH3 Monitoring




Ammonia calibration curve obtained using FTIR monitoring at 966cm<sup>-1</sup>. Good linearity is observed down to the low ppm range. The short detection pathlength given by the flow-tube geometry (16cm), however, makes it difficult to measure below 1ppm with confidence.




An improved detection limit for our apparatus can be achieved by converting ammonia into NO on a high temperature catalyst, and measuring the resulting change in background NO<sub>x</sub> concentration using a conventional NO<sub>x</sub> analyzer. This will allow us to work at more relevant ammonia concentrations.

## Observations



Reference spectrum of ammonium sulphate aerosol suspension. Total particle concentration:  $2 \times 10^5 \# \text{cm}^{-3}$  Mode: 150nm




Reference spectrum of ammonia. Ammonia partial pressure [NH<sub>3</sub>]: 20ppm.

## Results & Outlook

- Preliminary validation experiments show that our system can monitor both gas- and condensed phase composition.
- Good linearity is observed in gas-phase concentration changes.
- Ancillary deliquescence/efflorescence experiments show good reproducibility in terms of aerosol number density and composition.

Next proposed experimental steps include:

- Kinetic measurement of NH<sub>3</sub> uptake rate against aqueous sea-salt aerosols of different ionic composition.
- Heterogeneous reaction of NH<sub>3</sub> on sulfuric acid aerosols.
- Spectroscopic investigation of ammonium nitrate nucleation.



### Acknowledgements

Enterprise Ireland  
EU Marie Curie Fellowship  
HEA-PRTL cycle II Programme